

Richie Brace Fitting Guide



1) Unlace the shoe fully and slip the brace into shoe.



2) Slide the brace back so that it meets the heel counter in the back of the shoe.



3) Open the sides of the brace and slide your foot into the shoe.



4) Grab the sides of the brace and pull the brace to the back of the shoe.



5) Kick your heel to make sure the brace is at the back of the shoe.



6) Align the sides of the brace to the sides of your leg, making sure the brace runs up the middle of each side of your leg.



7) Secure the two front straps first.



8) Grab the back strap and feed it through the slot in the back of the brace. Then secure the back strap.



9) Push the brace back to the heel of the shoe one more time.



10) Tie the shoelaces.