



Account # _____

Name: _____
 Address: _____
 City: _____ State: _____ Zip _____
 Phone: _____ Email: _____

Freestyle OTC Orthotics

	A	B	C	D	E	F	G	H	I
Men's Sizes	4-4.5	5-5.5	6-6.5	7-7.5	8-8.5	9-9.5	10-11	12-13	14-15
Women's Sizes	6-6.5	7-7.5	8-8.5	9-9.5	10-10.5	11-11.5	12-13	14+	NA

<i>Freestyle Walker</i> (see back for product description)		Quantities									Total Units
4 Degrees Varus Ext RF Post	Full Length										
	Sulcus Length										
No RF Post	Full Length										
	Sulcus Length										
<i>Freestyle High Heel</i> (see back for product description)											
No RF Post	Full Length										
	Sulcus Length										
<i>Freestyle Support</i> (see back for product description)											
4 Degrees Varus Ext RF Post	Met Length										
	No Top Cover										
No RF Post	Met Length										
	No Top Cover										
<i>Freestyle Casual</i> (see back for product description)											
4 Degrees Varus Ext RF Post	Full Length										
	Sulcus Length										
No RF Post	Full Length										
	Sulcus Length										
<i>Freestyle Starter Kit</i> (see back for product description)											
Sizes A-E & Sample Kit	1 starter kit per practice/location	See back side for more details									



6299 Guion Rd., Indianapolis, IN 46268
 p. 800.444.3632 f. 800.233.2280 hours: M-Th: 7 a.m. - 5:30 p.m. EST
 e. customerservice@alliedosilabs.com w: www.alliedosilabs.com

Subtotal	
Shipping - 8%	
Total	



Freestyle OTC Orthotics

PRODUCT DESCRIPTIONS

<i>Freestyle Walker</i> (1/8" Triton Shell for Semiflexible Support)		<i>Freestyle High Heel</i> (1/8" Triton Shell for Semiflexible Support)		<i>Freestyle Support</i> (1/8" Triton Shell for Semiflexible Support)		<i>Freestyle Casual</i> (1/8" Triton Shell for Semiflexible Support)	
Use: Sports and Everyday Use		Use: Plantar Fasciitis & Heel Pain		Use: For Everyday Use		Use: Men's & Women's Dress Shoes	
Full or Sulcus Top Cover 1/8" EVA Padding	Non-Posted Rearfoot Optional	Full or Sulcus Top Cover 1/8" EVA Padding	No Shell In Heel - Hour Glass Shape For Better Shoe Fit	Met Top Cover or No Top Cover	Non-Posted Rearfoot Optional	Full or Sulcus Top Cover 1/8" EVA Padding	Hour Glass Shape & Hole in Heel For Better Shoe Fit
4 Degrees Varus Extrinsic Rearfoot Post	Suede Bottom Cover Under Forefoot Only	No Rearfoot Post	Suede Bottom Cover	4 Degrees Varus Extrinsic Rearfoot Post	No Suede Bottom Cover	4 Degrees Varus Extrinsic Rearfoot Post Non-Posted Optional	Suede Bottom Cover Under Forefoot Only

STANDARD Freestyle Starter Kit (20 Most Popular Sizes/Styles): \$375.00 plus FREE Samples

	A	B	C	D	E	F	G	H	I
Men's Sizes	4-4.5	5-5.5	6-6.5	7-7.5	8-8.5	9-9.5	10-11	12-13	14-15
Women's Sizes	6-6.5	7-7.5	8-8.5	9-9.5	10-10.5	11-11.5	12-13	14+	NA
	(1) Walker Full Length Posted	(2) Walker Full Length Posted	(2) Walker Full Length Posted	(2) Walker Full Length Posted	(2) Walker Full Length Posted	0	0	0	0
	(1) Support Met Length Posted	(3) Support Met Length Posted	(3) Support Met Length Posted	(3) Support Met Length Posted	(1) Support Met Length Posted	0	0	0	0

All Starter Kits Include a FREE Sample of Each Style

Freestyle Walker: Full-Length | Posted
Freestyle Casual: Sulcus-Length | No Post
Freestyle Support: Met-Length | Posted
Freestyle High Heel: Sulcus-Length | No Post

Contact Us And Ask How You Can BUILD YOUR OWN Starter Kit!



6299 Guion Rd., Indianapolis, IN 46268

p. 800.444.3632

f. 800.233.2280

hours: M-Th: 7: a.m. - 5:30 p.m. EST

e. customerservice@alliedosilabs.com

w: www.alliedosilabs.com